


High resolution, fast and more informative LA imaging with the icpTOF

O. Borovinskaya¹, M. Tanner¹, R. Belissont², M.C. Boiron², B. Luais², L. Summerfield³, S. Van Malderen⁴, T. Van Acker⁴, F. Vanhaecke⁴


¹TOFWERK, ²Université de Lorraine, ³ESI, ⁴University of Ghent
email: borovinskaya@tofwerk.com

- Fast and high resolution imaging is of growing interest for many geological and biological applications, but poses challenges to laser ablation (LA) units and mass spectrometers (MS).
- Recently developed Dual Concentric Injector (DCI) from ESI [1] and Aerosol Rapid Introduction System (Aris) from Teledyne [2] as add-ons to the existing ablation chambers enable rapid aerosol washout of <50 ms and increase imaging speed and signal to noise ratio.
- To utilize these new systems for multi-element imaging a mass spectrometer with high acquisition speed and simultaneous detection is required.
- This work demonstrates the coupling of TOFWERK time-of-flight mass spectrometer icpTOF with NRW 213 nm laser/TwoVol2 cell/DCI (ESI) and with Analyte G2 laser/Helix cell/Aris (Teledyne) for fast and easy imaging at high spatial resolution.


Imaging concept



Sphalerite mapping with NRW 213/TwoVol2 cell/DCI



Ceriodaphnia imaging with Analyte G2/Helix cell/Aris



Summary

- Rapid aerosol transfer systems resolve individual laser pulses at frequencies > 10 Hz, enabling sample scan at much higher speed than was earlier possible with standard configurations.
- icpTOF and laser synchronization on a single pixel basis simplify image reconstruction and provides sharp images with no pixel intermix.
- icpTOF measures several full mass spectra per every single pixel, providing higher information density about the sample.

Authors would like to acknowledge ESI and Teledyne for providing their laser ablation systems and for fruitful collaboration.

[1] D. N. Douglas et.al. Anal. Chem., 2015, 87, 11285-11294
[2] R. Belissont et.al. Geochimica et Cosmochimica Acta, 2014, 126, 518-540