

# TEKSCAN PRESSURE MEASUREMENT ACCURACY FOR ORTHOPAEDIC BIOMECHANICAL JOINT CONTACT MEASUREMENTS

Herregodts, S., ir.<sup>1,2</sup> – Verstraete, M.A., PhD<sup>1,2</sup> – Victor, J., MD, PhD<sup>1</sup>  
– De Baets, P., PhD<sup>2</sup>


Dept. of Orthopaedics and Traumatology, Ghent University

Dept. of Mechanical Construction and Production, Ghent University




**Multi-axial test rig properties**  
Compression testing  
Biaxial testing  
Adjustable degrees of freedom of specimens  
Knee implant mounting

**Multi-axial test rig capabilities**  
Research of the Tekscan sensor characteristics  
Tekscan sensor automated calibration  
Evaluation of knee implants: pressure distribution measurement and stability research




**Aim of this study:** Improve the accuracy and reliability of Tekscan intra-articular contact pressure measurements


## Sensor accuracy optimization technique



Automated sensor preconditioning and characterization  
Preconditioning cycle improves the sensor output stability  
Characterization of the sensor's step response allows for sensor drift compensating post-processing




**Automated multipoint calibration**  
Optimal contacting surfaces  
Fully automated and synchronized process



**Contact pressure measurements**



**Automated measurement data post-processing**  
Based on sensor specific behaviour and calibration data, sensor drift is eliminated during post-processing




Drift compensation technique reduces average error during realistic loading from 15,6 % to 3,4 %

## Sensor position measurement method




**Sensor fixation to the implant**





**Measurement method**

Implant specific 3D printed calibration block



**Sensor position computation**

Automated sensor position computation and visualization



**Pressure distribution visualisation**

Computation and visualization centres of pressure



Animated and synchronized visualization

Easy interpretation

## Conclusions

Optimized sensor preconditioning, sensor specific characterization and calibration combined with advanced data post-processing results in a 5 times higher measurement accuracy.

Sensor position measurement allows for intuitive visualization and interpretation of measurement data.



**Corresponding Author**  
Herregodts Stijn  
Ghent University  
Stijn.Herregodts@UGent.be